



# Homecare Product Catalogue

Laryngectomy and Tracheostomy Supplies



# Who Is NWOS Homecare?

## About Us

Homecare is an independent Dispensing Appliance Contractor (DAC), offering a wide range of products from all manufacturers.

We provide a home delivery service for stoma - ostomy, continence, tracheostomy and laryngectomy (head and neck) appliances.

We are completely independent, and for this reason, hold no vested interests in promoting or marketing any one brand of products directly to patients.

## Product Range

We stock a comprehensive range of products from all manufacturers in the following areas:



**Stoma - Ostomy**



**Head and Neck**



**Continence**



**Monika Menefel, Team Leader**



**Charlotte Atkinson, Assistant**



**Sharon Fisher, Assistant**

The Professional Services Team  
at NWOS Homecare

Across our ranges we work with a variety of manufacturers including:

- 3M
- Atos
- B Braun
- Bard
- Clini Supplies
- Coloplast
- Convatec
- CUI
- Dansac
- Flexicare
- Hollister
- Insight
- Oakmed
- Optimum
- Peak
- Pelican
- Salts
- Severn Healthcare
- Sheffmed
- SuportX
- Teleflex



## Why order from NWOS Homecare?

- + Independent company
- + Free delivery service
- + NHS electronic prescription

If you can't find what you need please call our friendly advisors who will be happy to help.



@ Email [homecareorders@nwossurgical.co.uk](mailto:homecareorders@nwossurgical.co.uk)

☎ Freephone **0800 243 103** 8.30am to 5.00pm

📱 SMS text number **07464 986 747**

www [www.homecarenwos.co.uk](http://www.homecarenwos.co.uk)



**Homecare**<sup>®</sup>  
from NWOS<sup>®</sup>

Baseplates	5
Heat and Moisture Exchange (HME) Cassettes	12
Tubes and Buttons	21
Valve Housing	29
Speaking Valves	30
Cleaning Devices	32
Dressings	38
Tube Holders	40
Protectors and Shields	43
Shower Protectors	47
Accessories	48



**Check with your Head and Neck Nurse, ENT Nurse or Speech Therapist for suitability**

**Please Note:** For all categories marked with this symbol, we recommend that you always check with your Head and Neck Nurse, ENT Nurse or Speech Therapist to ensure the best product for your requirements.

### Fahl Laryvox Tape Standard

The standard baseplate range offers excellent adhesive strength and are, transparent, water-repellent, available in 3 shapes.



Check with your Head and Neck Nurse, ENT Nurse or Speech Therapist for suitability



48100



48200

### Fahl Laryvox Tape Flexible

A unique and flexible base plate range with the highest adhesive strength, ideal for those individuals demonstrating issues with base plate retention or considering use of a hands free system. In addition to its exceptional level of adhesion the flexible base plate range are also; transparent, water-repellent, available in 3 sizes.



48220

	Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/>	48100	Fahl Laryvox Tape Standard	Round	15
<input type="checkbox"/>	48200	Fahl Laryvox Tape Standard	Oval	15
<input type="checkbox"/>	48300	Fahl Laryvox Tape Standard, XL	Oval	15
<input type="checkbox"/>	48120	Fahl Laryvox Tape Flexible	Round	15
<input type="checkbox"/>	48220	Fahl Laryvox Tape Flexible,	Oval	15
<input type="checkbox"/>	48320	Fahl Laryvox Tape Flexible, XL	Oval	15



### Fahl Laryvox Tape Hydrosoft

The Hydrosoft baseplate is made of soft hydrocolloid material and incorporates a very mild adhesive. Particularly skin friendly and comfortable on application. NOT SUITABLE FOR USE IN COMBINATION WITH SHOWER GUARD.

### Fahl Laryvox Tape Hypoallergen

A unique comfortable and flexible baseplate that has been designed for patients with skin sensitivities or skin reactions but also offers strong adhesion.

### Fahl Laryvox Tape Comfort

A comfortable tinted baseplate that adapts well to the contours of the skin and offers the strongest level of adhesion, making it the ideal choice for patients who have experienced issues with baseplate retention or wish to trial the hands free system e.g. the Blom-Singer ATSV II.

### Fahl Laryvox Tape Extra Fine

A unique transparent, strong, waterproof baseplate that provides long lasting adhesion. The ultra thin, transparent material used in the extra fine baseplate adapts well to all tracheostomal conditions, making it an ideal choice for those patients with an irregular or recessed stoma and peristomal areas.



48050



48130



48140

	Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/>	48130	Fahl Laryvox Tape Hydrosoft	Round	15
<input type="checkbox"/>	48230	Fahl Laryvox Tape Hydrosoft	Oval	15
<input type="checkbox"/>	48330	Fahl Laryvox Tape Hydrosoft	X Large oval	15
<input type="checkbox"/>	48050	Fahl Laryvox Tape Hypoallergen	Round	15
<input type="checkbox"/>	48150	Fahl Laryvox Tape Hypoallergen	Oval	15
<input type="checkbox"/>	48250	Fahl Laryvox Tape Hypoallergen	X Large oval	15
<input type="checkbox"/>	48140	Fahl Laryvox Tape Comfort	Round	15
<input type="checkbox"/>	48240	Fahl Laryvox Tape Comfort	Oval	15
<input type="checkbox"/>	48340	Fahl Laryvox Tape Comfort	X Large oval	15
<input type="checkbox"/>	48400	Fahl Laryvox Tape Extra Fine	-	15

### Blom-Singer TruSeal Contour Low Profile Adhesive Housing

It has a lower profile from the standard housing, reducing its height and overall profile.



BE6074



BE6075



BE6076

### Blom-Singer AccuFit Adhesive Housing

The Blom-Singer AccuFit® Adhesive Housing is designed to be the strongest adhesive option while still maintaining flexibility. It has increased elasticity and is less prone to tearing.

### HydroFit Hydrocolloid Adhesive Housing

The Blom-Singer HydroFit® Adhesive Housing is a hydrocolloid adhesive baseplate made with a protective gel-like layer and a moisture-absorbent adhesive. The HydroFit is gentle on the skin and flexible, but maintains strong adhesion. Available in three shapes: round, oval, and oval extra.

### Blom-Singer Regular Tape Discs


A variety of double-sided adhesive foam and tape discs to provide the attachment of the tracheostoma valve housing to the stoma.

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> BE6074	Blom-Singer TruSeal Contour Low Profile Adhesive Housing	Round	30
<input type="checkbox"/> BE6075	Blom-Singer TruSeal Contour Low Profile Adhesive Housing	Oval	30
<input type="checkbox"/> BE6076	Blom-Singer TruSeal Contour Low Profile Adhesive Housing	Standard	30
<input type="checkbox"/> BE6082	Blom-Singer AccuFit Adhesive Housing	Round	30
<input type="checkbox"/> BE6083	Blom-Singer AccuFit Adhesive Housing	Oval	30
<input type="checkbox"/> BE6084	Blom-Singer AccuFit Adhesive Housing	Oval Extra	30
<input type="checkbox"/> BE6085	Blom-Singer HydroFit Hydrocolloid Adhesive Housing	Round	30
<input type="checkbox"/> BE6086	Blom-Singer HydroFit Hydrocolloid Adhesive Housing	Oval	30
<input type="checkbox"/> BE6087	Blom-Singer HydroFit Hydrocolloid Adhesive Housing	Oval Extra	30
<input type="checkbox"/> BE6041	Blom-Singer Regular Tape Discs	Standard	30
<input type="checkbox"/> BE6042	Blom-Singer Regular Tape Discs	Large	30
<input type="checkbox"/> BE6043	Blom-Singer Thin Tape Discs	Standard	30
<input type="checkbox"/> BE6044	Blom-Singer Thin Tape Discs	Large	30
<input type="checkbox"/> BE6034	Blom-Singer Heavy Duty Tape Discs	Standard	30
<input type="checkbox"/> BE6035	Blom-Singer Heavy Duty Tape Discs	Large	30



### Provox Life Night Adhesive

Provox Life Night Adhesive is a skin-friendly hydrogel adhesive that soothes and cools the skin.



8261

### Provox Life Stability Adhesive

Provox Life Stability is a firm everyday adhesive that provides solidity, especially if your stoma is deep. It can be used if you plan on speaking a lot or if you like to speak hands-free as it gives you better stability with its specially designed base. It is also suitable for people with a deeper stoma.

### Provox Life Sensitive Adhesive

Provox Life Sensitive is a gentle everyday adhesive for sensitive skin. A comfortable, gentle and low-profile adhesive that forms a good seal and provides a secure fit.



**Check with your Head and Neck Nurse, ENT Nurse or Speech Therapist for suitability**

### Provox Life Standard Adhesive

Provox Life Standard is an adhesive for everyday use. It is soft, flexible and has a low-profile to provide a strong seal and secure fit around your stoma.



7466

#### **\* IMPORTANT INFORMATION:**

All Provox Life products can only be used with other Provox Life products as an integrated system. They cannot be used in conjunction with any similar products from other branded systems.

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> 8261	*Provox Life Night Adhesive	-	15
<input type="checkbox"/> 7463	*Provox Life Sensitive Adhesive	Round	30
<input type="checkbox"/> 7464	*Provox Life Sensitive Adhesive	Oval	30
<input type="checkbox"/> 7466	*Provox Life Sensitive Adhesive	Plus	30
<input type="checkbox"/> 7460	*Provox Life Standard Adhesive	Round	30
<input type="checkbox"/> 7461	*Provox Life Standard Adhesive	Oval	30
<input type="checkbox"/> 7462	*Provox Life Standard Adhesive	Plus	30
<input type="checkbox"/> 8263	*Provox Life Stability Adhesive	-	15

## Provox Flexiderm Adhesive

Offering strong, flexible, transparent adhesive properties.



7253

## Provox Optiderm Adhesive

Lighter adhesive properties designed to be very gentle on the skin.



7255

## Provox Xtrabase

A unique, convex self-adhering base plate that has a special shape and reinforced centre, ensuring it sits close to the stomal rim offering excellent retention. The XtraBase is ideal for people seeking strong adhesion and those with a deep/irregular shaped stoma.

## Provox StabiliBase

StabiliBase is a strong adhesive that provides better stability, especially if the stoma is deep. With new stabilising bars.

	Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/>	7253	Provox Flexiderm Adhesive	Round	20
<input type="checkbox"/>	7254	Provox Flexiderm Adhesive	Oval	20
<input type="checkbox"/>	7255	Provox Optiderm Adhesive	Round	20
<input type="checkbox"/>	7256	Provox Optiderm Adhesive	Oval	20
<input type="checkbox"/>	7265	Provox Xtrabase	Convex	20
<input type="checkbox"/>	7289	Provox StabiliBase	Convex	15
<input type="checkbox"/>	7318	Provox StabiliBase Optiderm	-	15
<input type="checkbox"/>	7331	Provox Flexiderm Adhesive	Plus	20
<input type="checkbox"/>	7332	Provox Optiderm Adhesive	Plus	20



### Spirotect Fabric Baseplate

Skin friendly fabric dressings that secure HME Cassettes at the laryngectomy stoma site. A large surface area for better adhesion, helping to prevent air leakage. Water resistant adhesive makes it suitable for showering. Can be used for up to 50 hours, meaning less changes and skin irritation.

### Spirotect Hydrocolloid Baseplate

Skin friendly, soft hydrocolloid material. Ideal for laryngectomees with sensitive skin.

### Spirotect Clear and Strong Baseplate

Discreet, transparent, lightweight and flexible design. Water resistant adhesive makes it suitable for showering.

PTW100



PTW200



PTW300



PTW101



PTW201



PTW301



Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> PTW100	Spirotect Fabric Baseplate	Oval	20
<input type="checkbox"/> PTW101	Spirotect Fabric Baseplate	Round	20
<input type="checkbox"/> PTW200	Spirotect Hydrocolloid Baseplate	Oval	20
<input type="checkbox"/> PTW201	Spirotect Hydrocolloid Baseplate	Round	20
<input type="checkbox"/> PTW300	Spirotect Clear & Strong Baseplate	Oval	20
<input type="checkbox"/> PTW301	Spirotect Clear & Strong Baseplate	Round	20
<input type="checkbox"/> PTW501	Spirotect Xtra Sticky Baseplate	Round	20
<input type="checkbox"/> PTW700	Spirotect Baseplate Secure Convex	-	20



**PHONIQ® DigiClose™ PhoniQ Baseplate**

PHONIQ® Baseplates adhere gently to the skin area around a tracheostoma to provide an airtight interface between the stoma and PHONIQ® DigiClose™ HME Filter. When used together inhaled air is effectively warmed, moistened and filtered to promote optimum lung health. PHONIQ® Baseplates can be worn both DAY and NIGHT for up to 3 days.

**Phoniq Baseplate Secure**

Made of Clear flexible material. Ensures effective seal for voicing and airway protection.



LAPBP0002

**Phoniq Baseplate Hydrocolloid**

Made of hydrocolloid material for sensitive skin. Thin and soft for improved contouring.



LAPBP0003

	Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/>	LAPBP0001	Phoniq Baseplate Secure	Oval	20
<input type="checkbox"/>	LAPBP0002	Phoniq Baseplate Secure	Round	20
<input type="checkbox"/>	LAPBP0003	Phoniq Baseplate Hydrocolloid	Oval	20
<input type="checkbox"/>	LAPBP0004	Phoniq Baseplate Hydrocolloid	Round	20



Fahl HME cassettes are discreet and low profile providing an effective filter to airborne particles entering the lungs.



**Check with your Head and Neck Nurse, ENT Nurse or Speech Therapist for suitability**

## Fahl Laryvox HME HighFlow

Low respiratory resistance therefore easy to breathe through, applicable for all patients following laryngectomy. Ideal starting point for all new users of HME's. Used during more strenuous physical exertion. Easy finger occlusion for speech.

## Fahl Laryvox O2 HME

Provides a unique combination of Laryvox HME with an O2 connection, ideal for the patient who is requiring supplementary oxygen.

## Fahl Laryvox HME

Normal respiratory resistance, applicable for all patients following laryngectomy. Used for normal physical exertion, easy finger occlusion for speech.



49800



49802



49810

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> 46460	Fahl Humidotwin HME	-	30
<input type="checkbox"/> 49800	Fahl Laryvox HME	-	30
<input type="checkbox"/> 49810	Fahl Laryvox HME HighFlow	-	30
<input type="checkbox"/> 49802	Fahl Laryvox O2 HME	-	30
<input type="checkbox"/> 48460	Fahl Laryvox Pro HME HighFlow cassette	-	30
<input type="checkbox"/> 48461	Fahl Laryvox PRO HME Cassette	-	30

## Fahl Laryvox Touch HME Highflow

A heat and moisture Exchange (HME) cassette range with low pressure finger occlusion, A highFlow HME range with low respiratory resistance, suitable for new users of HME cassettes or during physical exertion.

49850-02



49850-04



49850-13



49862



## Fahl Laryvox Extra HME


The first HME to be available in 3 different resistance levels, it filters, warms and humidifies respiratory air. An easy to operate speaking option for patients with voice prosthesis, with an unobstructive, inconspicuous design.

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> 49850-01	Fahl Laryvox Touch, Blue Metallic	-	30
<input type="checkbox"/> 49850-02	Fahl Laryvox Touch, Gold	-	30
<input type="checkbox"/> 49850-03	Fahl Laryvox Touch, Silver	-	30
<input type="checkbox"/> 49850-04	Fahl Laryvox Touch, Beige	-	30
<input type="checkbox"/> 49850-13	Fahl Laryvox Touch, Clear,	-	30
<input type="checkbox"/> 49860	Fahl Laryvox Extra HME HighFlow	-	30
<input type="checkbox"/> 49862	Fahl Laryvox Extra HME Medium	-	30
<input type="checkbox"/> 49861	Fahl Laryvox Extra HME HighFlow	-	30
<input type="checkbox"/> 49863	Fahl Laryvox Extra HME Sport	-	30



## Provox Life Go HME

Provox Life Go HME provides lower breathing resistance. Provox Life Go HME is designed to make it easier to breathe when you are out and about or engaging in daily activities. This HME has a low profile and a smooth rounded form for a more discreet look.



## Provox Life Home HME

Provox Life Home HME has the highest humidification level of the daytime HMEs. Moreover, Provox Life Home HME is designed for at-home activities, such as reading.

## Provox Life Energy HME

Provox Life Energy HME provides good humidification and low breathing resistance. Designed to support you when you are being physically active, such as running or cycling.

## Provox Life Night HME

Provox Life Night HME helps reduce coughing for more restful nights. Provox Life Night HME offers high humidification which helps reduce mucus production and coughing, allowing you to rest while improving your lung health.

## Provox Life Protect HME

Provox Life Protect HME is designed for use when you are away from home and need extra protection, such as travelling on buses or trains; or when shopping or spending time outdoors.

**\* IMPORTANT INFORMATION:** All Provox Life products can only be used with other Provox Life products as an integrated system. They cannot be used in conjunction with any similar products from other branded systems.

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> 7440	*Provox Life Freehands HME	-	30
<input type="checkbox"/> 8310	*Provox Life Go HME cassettes	-	30
<input type="checkbox"/> 8311	*Provox Life Home HME cassettes	-	30
<input type="checkbox"/> 8312	*Provox Life Energy HME cassettes	-	30
<input type="checkbox"/> 8262	*Provox Life Night HME cassettes	-	30
<input type="checkbox"/> 8313	*Provox Life Protect HME cassettes	-	15

**Provox Micron HME**

A highly effective bacterial and viral filter (over 99%). Ideal for use in crowded places such as schools, supermarkets and on public transport.



**Check with your Head and Neck Nurse, ENT Nurse or Speech Therapist for suitability**



7248

**Provox Xtramoist HME**

Helps to keep lungs healthy by humidifying and filtering the air the user breathes. To be used at night to get the best possible humidification and helping the user breathe comfortably.

**Provox Xtraflow**

Helps to keep lungs healthy by humidifying and filtering the air the user breathes. To be used during the day when the user is physically active, as it has the lowest breathing resistance.



7290



7291

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> 7248	Provox Micron HME	-	30
<input type="checkbox"/> 7290	Provox Xtramoist HME	-	30
<input type="checkbox"/> 7291	Provox Xtraflow	-	30

PHONIQ® DigiClose™ HME Filters feature a patented pressure release mechanism which activates during reflexive coughing or sneezing. PHONIQ® DigiClose™ Active and Comfort are available in two colours beige and clear.



**Check with your Head and Neck Nurse, ENT Nurse or Speech Therapist for suitability**

## PHONIQ® DigiClose™ Active Filter

The PHONIQ® DigiClose™ Active provides low breathing resistance and allows for high flow rates for exercise or possibly those with lung disease.



LAPFT0002

## PHONIQ® Low Profile HME Filter

The PHONIQ® Low Profile HME fits perfectly in the PHONIQ® stoma buttons, laryngectomy tubes, and baseplates available. Discreet low-profile design. Allows regular flow and resistance for normal everyday activities.



LAPFT0005

## PHONIQ® DigiClose™ Comfort Filter

The PHONIQ® DigiClose™ Comfort allows regular flow and resistance for normal every day activity levels.

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> LAPFT0002	PhonIQ DigiClose Active HME Filter (Clear)	-	30
<input type="checkbox"/> LAPFT0004	PhonIQ DigiClose Active HME Filter (Beige)	-	30
<input type="checkbox"/> LAPFT0005	PhonIQ Low Profile HME Filter (Clear)	-	30
<input type="checkbox"/> LAPFT0001	PhonIQ DigiClose Comfort HME Filter (Clear)	-	30
<input type="checkbox"/> LAPFT0003	PhonIQ DigiClose Comfort HME Filter (Beige)	-	30

### Spirotect Hands Free

Antibacterial coated foam to help reduce the potential of respiratory lung infection. Hands free HME enabling the laryngectomee to speak without using fingers.



**Check with your Head and Neck Nurse, ENT Nurse or Speech Therapist for suitability**



PTC1004

#### Normal

Suitable for day or night use. It allows the laryngectomee to easily obtain airflow closure and speech by using a light finger touch.



PTC1005

#### Lateral

Suitable for day or night use. It allows the laryngectomee to easily obtain airflow closure and speech by gently pushing the top of the cassette. The size of the finger will not affect the closure effectiveness.

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> PTC1004	Spirotect Hands Free Anti Bacterial Cassettes Normal	-	30
<input type="checkbox"/> PTC1005	Spirotect Hands Free Anti Bacterial Cassettes Lateral	-	30

### Blom-Singer MHE Cartridge Easy Flow

With EasyFlow, you can breathe a little easier. Its foam filter is more porous than the ClassicFlow. This allows you to breathe more freely to match your activity level and pulmonary needs. The beige tint and non-reflective surface texture make the EasyFlow less noticeable.



**Check with your Head and Neck Nurse, ENT Nurse or Speech Therapist for suitability**



### Blom-Singer HumidiFilter Holder

The HumidiFilter Holder is a lightweight device that attaches to the front of a tracheostoma adhesive housing.

### Blom-Singer HumidiFilter Holder Foam Filter

For use with Blom-Singer HumidiFilter Holder (BE1060).

### Blom-Singer Speakfree HME EasyFlow

With EasyFlow the foam filter is more porous than the ClassicFlow. This allows patient to breath more freely to match their activity level and pulmonary needs.

### Blom-Singer Speakfree HME ClassicFlow

ClassicFlow helps assure you breathe warmed, filtered and humidified air. It is designed for increased moisture retention and resistance.

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> BE1055EF	Blom-Singer HME Cartridge EasyFlow	-	30
<input type="checkbox"/> BE1055EZ	Blom-Singer HME Cartridge Classic Flow	-	30
<input type="checkbox"/> BE1060	Blom-Singer HumidiFilter Holder	-	1
<input type="checkbox"/> BE1080	Blom-Singer HumidiFilter Holder Foam Filter	-	30
<input type="checkbox"/> BE1090EF	Blom-Singer Speakfree HME EasyFlow	-	30
<input type="checkbox"/> BE1090EZ	Blom-Singer Speakfree HME ClassicFlow	-	30
<input type="checkbox"/> BE1091	Blom-Singer Day&Night HME Cartridge ClassicFlow	-	30
<input type="checkbox"/> BE1092	Blom-Singer Day&Night HME Cartridge EasyFlow	-	30

**Freevent DualCare HME**

Freevent DualCare combines a speaking valve with a HME enabling speaking and supporting your lungs with important heating and moisturising.



7742

 Check with your Head and Neck Nurse, ENT Nurse or Speech Therapist for suitability

**Provox Trachphone**

A HME offering effective pulmonary and vocal rehabilitation for tracheostomy patients.



7707

7767/8



**Freevent XtraCare HME**

Freevent XtraCare combines an HME with a filter that provides more than humidification. It also provides very effective filtration, for extra protection.

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> 7742	Freevent DualCare HME	-	30
<input type="checkbox"/> 7707	Trachphone	-	30
<input type="checkbox"/> 7767	Freevent XtraCare HME (white)	-	30
<input type="checkbox"/> 7768	Freevent XtraCare HME (blue)	-	30

## Provox Freehands FlexiVoice

Offers both the freedom to speak hands-free and with finger occlusion when wanted. For hands-free speech the user can choose between three different strengths: light, medium, and strong. The choice depends on the amount of airflow and pressure required for speaking and breathing.



8013



8161

## Provox Filter Cassette

A disposable HME cassette that simply clicks onto the base of the Provox Freehands System, offering consistent conditioning of inhaled air. The HME also protects the Freehands device from secretions, whilst providing an effective level of airflow resistance.

Product Code	Item Description	Size/Shape	Pack Size	
<input type="checkbox"/>	8161	Provox Freehands FlexiVoice Light	-	1
<input type="checkbox"/>	8162	Provox Freehands FlexiVoice Medium	-	1
<input type="checkbox"/>	8163	Provox Freehands FlexiVoice Strong	-	1
<input type="checkbox"/>	8221	Provox FreeHands Flow	-	30
<input type="checkbox"/>	8220	Provox FreeHands Moist	-	30
<input type="checkbox"/>	8013	Provox Luna HME	-	30
<input type="checkbox"/>	8014	Provox Luna Adhesive	-	15
<input type="checkbox"/>	8021	Provox Freehands Support	Flat	1
<input type="checkbox"/>	8022	Provox Freehands Support	Medium	1
<input type="checkbox"/>	8023	Provox Freehands Support	Deep	1
<input type="checkbox"/>	8024	Provox Freehands Support Adhesive	-	15

## Phoniq Laryngectomy Tube

Transparent high-quality medical grade silicone tube, universal 22mm connection. Available in fenestrated and non-fenestrated option.



Check with your Head and Neck Nurse, ENT Nurse or Speech Therapist for suitability

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> LAPLF0840-1440	Phoniq Laryngectomy Tube (Fenestrated) 40mm length	Size 8-14	1
<input type="checkbox"/> LAPLF0850-1450	Phoniq Laryngectomy Tube (Fenestrated) 50mm length	Size 8-14	1
<input type="checkbox"/> LAPLF0860-1460	Phoniq Laryngectomy Tube (Fenestrated) 60mm length	Size 8-14	1

## Blom-Singer StomaSoft Laryngectomy Tube

Blom-Singer® StomaSoft® Laryngectomy Tubes are flexible, soft, and compatible with HME and other accessories from InHealth Technologies®. They feature reinforced strap holders and are available with or without fenestrations

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> BE6698	Blom-Singer StomaSoft Laryngectomy Tube (Non - Fenestrated)	8/36	1
<input type="checkbox"/> BE6699	Blom-Singer StomaSoft Laryngectomy Tube (Non - Fenestrated)	8/55	1
<input type="checkbox"/> BE6702	Blom-Singer StomaSoft Laryngectomy Tube (Non - Fenestrated)	10/36	1
<input type="checkbox"/> BE6703	Blom-Singer StomaSoft Laryngectomy Tube (Non - Fenestrated)	10/55	1
<input type="checkbox"/> BE6698F	Blom-Singer StomaSoft Laryngectomy Tube (Fenestrated)	8/36	1
<input type="checkbox"/> BE6699F	Blom-Singer StomaSoft Laryngectomy Tube (Fenestrated)	8/55	1
<input type="checkbox"/> BE6702F	Blom-Singer StomaSoft Laryngectomy Tube (Fenestrated)	10/36	1
<input type="checkbox"/> BE6703F	Blom-Singer StomaSoft Laryngectomy Tube (Fenestrated)	10/55	1



### Provox Life LaryTube Fenestrated

The Provox LaryTube can be used to attach a HME or a Provox FreeHands Speaking Valve. Provox LaryTube is available in Standard, Fenestrated and With Ring (to be secured with Provox adhesive).



7434

**\* IMPORTANT INFORMATION:**  
 All Provox Life products can only be used with other Provox Life products as an integrated system. They cannot be used in conjunction with any similar products from other branded systems.

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> 7429	*Provox Life LaryTube (Fenestrated) 36mm	8	1
<input type="checkbox"/> 7431	*Provox Life LaryTube (Fenestrated) 36mm	9	1
<input type="checkbox"/> 7433	*Provox Life LaryTube (Fenestrated) 36mm	10	1
<input type="checkbox"/> 7435	*Provox Life LaryTube (Fenestrated) 36mm	12	1
<input type="checkbox"/> 7430	*Provox Life LaryTube (Fenestrated) 55mm	8	1
<input type="checkbox"/> 7432	*Provox Life LaryTube (Fenestrated) 55mm	9	1
<input type="checkbox"/> 7434	*Provox Life LaryTube (Fenestrated) 55mm	10	1
<input type="checkbox"/> 7436	*Provox Life LaryTube (Fenestrated) 55mm	12	1
<input type="checkbox"/> 8048	*Provox Life LaryTube (Fenestrated) with Ring 36mm	8	1
<input type="checkbox"/> 8050	*Provox Life LaryTube (Fenestrated) with Ring 36mm	9	1
<input type="checkbox"/> 8052	*Provox Life LaryTube (Fenestrated) with Ring 36mm	10	1
<input type="checkbox"/> 8054	*Provox Life LaryTube (Fenestrated) with Ring 36mm	12	1
<input type="checkbox"/> 8049	*Provox Life LaryTube (Fenestrated) with Ring 55mm	8	1
<input type="checkbox"/> 8051	*Provox Life LaryTube (Fenestrated) with Ring 55mm	9	1
<input type="checkbox"/> 8053	*Provox Life LaryTube (Fenestrated) with Ring 55mm	10	1
<input type="checkbox"/> 8055	*Provox Life LaryTube (Fenestrated) with Ring 55mm	12	1



### Provox Life LaryTube Standard with Ring

The Provox LaryTube can be used to attach a HME or a Provox FreeHands Speaking Valve. Provox LaryTube is available in Standard, Fenestrated and With Ring (to be secured with Provox adhesive).



7422

**\* IMPORTANT INFORMATION:**

All Provox Life products can only be used with other Provox Life products as an integrated system. They cannot be used in conjunction with any similar products from other branded systems.

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> 7409	*Provox Life LaryTube Standard 27mm	8	1
<input type="checkbox"/> 7412	*Provox Life LaryTube Standard 27mm	9	1
<input type="checkbox"/> 7415	*Provox Life LaryTube Standard 27mm	10	1
<input type="checkbox"/> 7418	*Provox Life LaryTube Standard 27mm	12	1
<input type="checkbox"/> 7410	*Provox Life LaryTube Standard 36mm	8	1
<input type="checkbox"/> 7413	*Provox Life LaryTube Standard 36mm	9	1
<input type="checkbox"/> 7416	*Provox Life LaryTube Standard 36mm	10	1
<input type="checkbox"/> 7419	*Provox Life LaryTube Standard 36mm	12	1
<input type="checkbox"/> 7411	*Provox Life LaryTube Standard 55mm	8	1
<input type="checkbox"/> 7414	*Provox Life LaryTube Standard 55mm	9	1
<input type="checkbox"/> 7417	*Provox Life LaryTube Standard 55mm	10	1
<input type="checkbox"/> 7420	*Provox Life LaryTube Standard 55mm	12	1
<input type="checkbox"/> 7421	*Provox Life LaryTube with Ring 36mm	8	1
<input type="checkbox"/> 7423	*Provox Life LaryTube with Ring 36mm	9	1
<input type="checkbox"/> 7425	*Provox Life LaryTube with Ring 36mm	10	1
<input type="checkbox"/> 7427	*Provox Life LaryTube with Ring 36mm	12	1
<input type="checkbox"/> 7422	*Provox Life LaryTube with Ring 55mm	8	1
<input type="checkbox"/> 7424	*Provox Life LaryTube with Ring 55mm	9	1
<input type="checkbox"/> 7426	*Provox Life LaryTube with Ring 55mm	10	1
<input type="checkbox"/> 7428	*Provox Life LaryTube with Ring 55mm	12	1



### Laryngotec Kombi Laryngectomy Tube

The Laryngotec Kombi Lingo silicone laryngectomy and tracheostomy tube is made of soft and flexible silicone for laryngectomees with a voice prosthesis and tracheostomees.

- Designed for long-term use, can be used up to 3-6 months.
- Medical grade silicone: soft, pliable, lightweight.
- Neck flange perfectly adapted to neck anatomy, with lateral openings for fixation of tube holder.
- Transparent, inconspicuous design.
- With 22 mm adapter (KOMBI).
- Each size in three different lengths.
- Custom-made sieve-fenestration possible.

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/>	14920-0744 Laryngotec Kombi Laryngectomy Tube (Non Fenestrated) 44mm	7	1
<input type="checkbox"/>	14920-0844 Laryngotec Kombi Laryngectomy Tube (Non Fenestrated) 44mm	8	1
<input type="checkbox"/>	14920-0944 Laryngotec Kombi Laryngectomy Tube (Non Fenestrated) 44mm	9	1
<input type="checkbox"/>	14920-1044 Laryngotec Kombi Laryngectomy Tube (Non Fenestrated) 44mm	10	1
<input type="checkbox"/>	14920-1144 Laryngotec Kombi Laryngectomy Tube (Non Fenestrated) 44mm	11	1
<input type="checkbox"/>	14920-1244 Laryngotec Kombi Laryngectomy Tube (Non Fenestrated) 44mm	12	1
<input type="checkbox"/>	14920-1344 Laryngotec Kombi Laryngectomy Tube (Non Fenestrated) 44mm	13	1
<input type="checkbox"/>	14920-0762 Laryngotec Kombi Laryngectomy Tube (Non Fenestrated) 62mm	7	1
<input type="checkbox"/>	14920-0862 Laryngotec Kombi Laryngectomy Tube (Non Fenestrated) 62mm	8	1
<input type="checkbox"/>	14920-0962 Laryngotec Kombi Laryngectomy Tube (Non Fenestrated) 62mm	9	1
<input type="checkbox"/>	14920-1062 Laryngotec Kombi Laryngectomy Tube (Non Fenestrated) 62mm	10	1
<input type="checkbox"/>	14920-1162 Laryngotec Kombi Laryngectomy Tube (Non Fenestrated) 62mm	11	1
<input type="checkbox"/>	14920-1262 Laryngotec Kombi Laryngectomy Tube (Non Fenestrated) 62mm	12	1
<input type="checkbox"/>	14920-1362 Laryngotec Kombi Laryngectomy Tube (Non Fenestrated) 62mm	13	1
<input type="checkbox"/>	14920-0765 Laryngotec Kombi Laryngectomy Tube (Non Fenestrated) 65mm	7	1
<input type="checkbox"/>	14920-0870 Laryngotec Kombi Laryngectomy Tube (Non Fenestrated) 70mm	8	1
<input type="checkbox"/>	14920-0975 Laryngotec Kombi Laryngectomy Tube (Non Fenestrated) 75mm	9	1
<input type="checkbox"/>	14920-1085 Laryngotec Kombi Laryngectomy Tube (Non Fenestrated) 85mm	10	1
<input type="checkbox"/>	14920-1187 Laryngotec Kombi Laryngectomy Tube (Non Fenestrated) 87mm	11	1
<input type="checkbox"/>	14920-1290 Laryngotec Kombi Laryngectomy Tube (Non Fenestrated) 90mm	12	1
<input type="checkbox"/>	14920-1390 Laryngotec Kombi Laryngectomy Tube (Non Fenestrated) 90mm	13	1



## Laryngotec Kombi Lingo Laryngectomy Tube

The Laryngotec Kombi Lingo silicone laryngectomy and tracheostomy tube is made of soft and flexible silicone for laryngectomees with a voice prosthesis and tracheostomees.

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/>	14720-0744 Laryngotec Kombi Lingo Laryngectomy Tube (Fenestrated) 44mm	7	1
<input type="checkbox"/>	14720-0844 Laryngotec Kombi Lingo Laryngectomy Tube (Fenestrated) 44mm	8	1
<input type="checkbox"/>	14720-0944 Laryngotec Kombi Lingo Laryngectomy Tube (Fenestrated) 44mm	9	1
<input type="checkbox"/>	14720-1044 Laryngotec Kombi Lingo Laryngectomy Tube (Fenestrated) 44mm	10	1
<input type="checkbox"/>	14720-1144 Laryngotec Kombi Lingo Laryngectomy Tube (Fenestrated) 44mm	11	1
<input type="checkbox"/>	14720-1244 Laryngotec Kombi Lingo Laryngectomy Tube (Fenestrated) 44mm	12	1
<input type="checkbox"/>	14720-1344 Laryngotec Kombi Lingo Laryngectomy Tube (Fenestrated) 44mm	13	1
<input type="checkbox"/>	14720-0762 Laryngotec Kombi Lingo Laryngectomy Tube (Fenestrated) 62mm	7	1
<input type="checkbox"/>	14720-0862 Laryngotec Kombi Lingo Laryngectomy Tube (Fenestrated) 62mm	8	1
<input type="checkbox"/>	14720-0962 Laryngotec Kombi Lingo Laryngectomy Tube (Fenestrated) 62mm	9	1
<input type="checkbox"/>	14720-1062 Laryngotec Kombi Lingo Laryngectomy Tube (Fenestrated) 62mm	10	1
<input type="checkbox"/>	14720-1162 Laryngotec Kombi Lingo Laryngectomy Tube (Fenestrated) 62mm	11	1
<input type="checkbox"/>	14720-1262 Laryngotec Kombi Lingo Laryngectomy Tube (Fenestrated) 62mm	12	1
<input type="checkbox"/>	14720-1362 Laryngotec Kombi Lingo Laryngectomy Tube (Fenestrated) 62mm	13	1
<input type="checkbox"/>	14720-0765 Laryngotec Kombi Lingo Laryngectomy Tube (Fenestrated) 65mm	7	1
<input type="checkbox"/>	14720-0870 Laryngotec Kombi Lingo Laryngectomy Tube (Fenestrated) 70mm	8	1
<input type="checkbox"/>	14720-0975 Laryngotec Kombi Lingo Laryngectomy Tube (Fenestrated) 75mm	9	1
<input type="checkbox"/>	14720-1085 Laryngotec Kombi Lingo Laryngectomy Tube (Fenestrated) 85mm	10	1
<input type="checkbox"/>	14720-1187 Laryngotec Kombi Lingo Laryngectomy Tube (Fenestrated) 87mm	11	1
<input type="checkbox"/>	14720-1290 Laryngotec Kombi Lingo Laryngectomy Tube (Fenestrated) 90mm	12	1
<input type="checkbox"/>	14720-1390 Laryngotec Kombi Lingo Laryngectomy Tube (Fenestrated) 90mm	13	1



## Provox LaryTube

The Provox LaryTube can be used to attach a HME or a Provox FreeHands Speaking Valve. Provox LaryTube is available in Standard, Fenestrated and With Ring (to be secured with Provox adhesive).



Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> 7637/7640/7643/7646	Provox Larytube (Fenestrated) 36mm length	Size 8-14	1
<input type="checkbox"/> 7638/7641/7644/7647	Provox Larytube (Fenestrated) 55mm length	Size 8-14	1
<input type="checkbox"/> 7601/7605/7609/7613	Provox Larytube Standard 27mm length	Size 8-14	1
<input type="checkbox"/> 7602/7606/7610/7614	Provox Larytube Standard 36mm length	Size 8-14	1
<input type="checkbox"/> 7603/7607/7611/7615	Provox Larytube Standard 55mm length	Size 8-14	1
<input type="checkbox"/> 7624/7626/7628/7630	Provox Larytube with Ring 36mm length	Size 8-14	1
<input type="checkbox"/> 7625/7627/7629/7631	Provox Larytube with Ring 55mm length	Size 8-14	1

## Duratwix Uni-Vario Lingo-Cuff Tracheostomy Tube

Tracheostomy tube with two inner cannulas, both with 15 mm connector, one with rotating 15 mm swivel connector (ICV), the other one fenestrated (IFCU).



Check with your Head and Neck Nurse, ENT Nurse or Speech Therapist for suitability

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> 19661	Duratwix Uni-Vario Lingo-Cuff tracheostomy tube with 2 inner cannulas - Standard	Size 6-10	1
<input type="checkbox"/> 19662	Duratwix Uni-Vario Lingo-Cuff tracheostomy tube with 2 inner cannulas - Short	Size 7-10	1
<input type="checkbox"/> 19211	Duratwix Uni-Vario tracheostomy tube with 2 inner cannulas- Standard	Size 6-10	1
<input type="checkbox"/> 19212	Duratwix Uni-Vario tracheostomy tube with 2 inner cannulas - Short	Size 7-10	1
<input type="checkbox"/> 19261	Duratwix Uni-Vario Lingo tracheostomy tube with 2 inner cannulas - Standard	Size 6-10	1
<input type="checkbox"/> 19262	Duratwix Uni-Vario Lingo tracheostomy tube with 2 inner cannulas - Short	Size 7-10	1
<input type="checkbox"/> 19611	Duratwix Uni-Vario Cuff tracheostomy tube with 2 inner cannulas - Standard	Size 6-10	1
<input type="checkbox"/> 19612	Duratwix Uni-Vario Cuff tracheostomy tube with 2 inner cannulas - Short	Size 7-10	1

### Provox Larybutton

The Provox LaryButton is a silicone stoma button, that maintains the opening of the stoma. It can be used to attach an Heat and Moisture Exchanger (HME), also in combination with the Provox FreeHands FlexiVoice speaking valve. The Provox LaryButton is well suited for patients with a circumferential skin fold around the stoma (a so-called “lip”) and patients with sensitive skin.



**Check with your Head and Neck Nurse, ENT Nurse or Speech Therapist for suitability**



Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> 7671	Provox Larybutton 8mm length, 12mm outer	12/8	1
<input type="checkbox"/> 7672	Provox Larybutton 8mm length, 14mm outer	14/8	1
<input type="checkbox"/> 7673	Provox Larybutton 8mm length, 16mm outer	16/8	1
<input type="checkbox"/> 7674	Provox Larybutton 8mm length, 18mm outer	18/8	1
<input type="checkbox"/> 7685	Provox Larybutton 18mm length, 12mm outer	12/18	1
<input type="checkbox"/> 7686	Provox Larybutton 18mm length, 14mm outer	14/18	1
<input type="checkbox"/> 7687	Provox Larybutton 18mm length, 16mm outer	16/18	1
<input type="checkbox"/> 7688	Provox Larybutton 18mm length, 18mm outer	18/18	1

### Laryngotec Kombi Stoma Buttons

The LARYNGOTEC® stoma buttons, made of soft, flexible silicone, are designed to help keep the tracheostoma open and can be used in conjunction with a Heat and Moisture Exchange (HME) cassette e.g the Blom-Singer EasyTouch HME System or a Hands-Free valve system such as the Blom-Singer ATSV II.

	Product Code	Item Description	Size/Shape	Pack Size
	15950-08	Laryngotec Kombi Stoma Button Fix 22mm	8	1
	15950-09	Laryngotec Kombi Stoma Button Fix 22mm	9	1
	15950-10	Laryngotec Kombi Stoma Button Fix 22mm	10	1
	15950-11	Laryngotec Kombi Stoma Button Fix 22mm	11	1
	15950-12	Laryngotec Kombi Stoma Button Fix 22mm	12	1
	15960-08	Laryngotec Kombi Stoma Button Fix 15mm	8	1
	15960-09	Laryngotec Kombi Stoma Button Fix 15mm	9	1
	15960-10	Laryngotec Kombi Stoma Button Fix 15mm	10	1
	15960-11	Laryngotec Kombi Stoma Button Fix 15mm	11	1
	15960-12	Laryngotec Kombi Stoma Button Fix 15mm	12	1

### Blom-Singer Tracheostoma Valve Housing

The Tracheostoma Valve Housing is a reusable, flexible housing that attaches the ATSV II or HME cartridge over the stoma.



Check with your Head and Neck Nurse, ENT Nurse or Speech Therapist for suitability



BE6038



BE6039

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> BE6038	Blom-Singer Tracheostoma Valve Housing	Standard	1
<input type="checkbox"/> BE6039	Blom-Singer Tracheostoma Valve Housing	Large	1

### Blom-Singer Adjustable Tracheostoma Valve II

The ATSV II consists of two pieces, a diaphragm/faceplate and body. The assembled two-piece valve inserts into the tracheostoma adhesive housing that is applied to the skin surrounding the stoma.



**Check with your Head and Neck Nurse, ENT Nurse or Speech Therapist for suitability**



BE8025H

### Blom-Singer Replacement Diaphragm/Faceplate

Replacement diaphragm/faceplate for Blom-Singer ATSV II Valve.

### Blom-Singer ATSV II Cap and 7 Foam Filters

The ATSV II HumidiFilter Cap is a twist-on plastic filter cap that fits over the ATSV II. It is designed to provide the laryngectomy with warmed, humidified and filtered air.

### Blom-Singer ATSV II HumidiFilter Cap Foam Filters

Replacement foam filters for use with ATSV II HumidiFilter Cap (BE1010).

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> BE8025H	Blom-Singer Adjustable Tracheostoma Valve II	-	1
<input type="checkbox"/> BE8026H	Blom-Singer Replacement Diaphragm/Faceplate	-	1
<input type="checkbox"/> BE1010	Blom-Singer ATSV II Cap and 7 Foam Filters	-	1
<input type="checkbox"/> BE1030	Blom-Singer ATSV II HumidiFilter Cap Foam Filters	-	30

## Fahl Tracheostomy Speaking Valve

The Tracheostomy Speaking Valve is developed to help with the treatment of Tracheostomy patients. The valve provides patients the ability for voicing and to communicate more effectively. Speaking valves should never be used with the tracheostomy tube cuff inflated as this can lead to asphyxiation and respiratory arrest.



Check with your Head and Neck Nurse, ENT Nurse or Speech Therapist for suitability



TSV/100

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> 27131	Fahl Combiphon Speaking Valve	-	1
<input type="checkbox"/> TSV/100	Fahl Trachie-Voice Tracheostomy Speaking Valve	-	2
<input type="checkbox"/> 7744	Freevent DualCare Speaking Valve (tinted)	-	1
<input type="checkbox"/> 7755	Freevent DualCare Speaking Valve (blue)	-	1
<input type="checkbox"/> TRVSV2000	Voiceline li Speaking Valve	White	2
<input type="checkbox"/> TRVSV2001	Voiceline li Speaking Valve	Green	2
<input type="checkbox"/> TRVSV2002	Voiceline li Speaking Valve	Orange	2
<input type="checkbox"/> TRVSV2003	Voiceline li Speaking Valve	Pink	2
<input type="checkbox"/> TRVSV2004	Voiceline li Speaking Valve	Blue	2




### Kapitex Tracheo Cleaning Brush

Kapitex cleaning brushes are used to ensure that the outer and inner tubes are kept free of dry mucus and debris. The brushes are not suitable for Tracheo PVC tubes or silicone tubes.



**Check with your Head and Neck Nurse, ENT Nurse or Speech Therapist for suitability**



Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> TRACC0001	Tracheo Cleaning Brush	8mm	2
<input type="checkbox"/> TRACC0002	Tracheo Cleaning Brush	10mm	2
<input type="checkbox"/> TRACC0003	Tracheo Cleaning Brush	12mm	2
<input type="checkbox"/> TRACC0004	Tracheo Cleaning Brush	14mm	2
<input type="checkbox"/> 7660	Provox TubeBrush	8mm	6
<input type="checkbox"/> 7661	Provox TubeBrush	12mm	6



### Laryvox Duo Brush

The Laryvox® Duo Brush is intended for daily cleaning of the voice prosthesis.

- Two in One cleaning system – sponge and length specific cleaning brush.
- Packaged in the SmartBox.
- External Mirror mounted to the SmartBox.
- Comes with a clip to attach a torch (sold separately).
- The length of the Laryvox® Duo Brush is chosen to correspond to different sizes of voice prostheses.
- Supplied in Laryvox® Duo Brush for voice prostheses in cases containing six pieces.



29110-08

	Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/>	29110-06	Laryvox Duo Brush	4-6mm	6
<input type="checkbox"/>	29110-08	Laryvox Duo Brush	7-8mm	6
<input type="checkbox"/>	29110-10	Laryvox Duo Brush	9-10mm	6
<input type="checkbox"/>	29110-12	Laryvox Duo Brush	11-13mm	6



## Blom Singer 16/20 Fr

An integral part of daily maintenance. The brush has important safety features to prevent it from being inserted too far into the prosthesis or trachea.



IN4064



Check with your Head and Neck Nurse, ENT Nurse or Speech Therapist for suitability

## Insight

Manufactured with soft nylon bristles that prevent damaging the silicone inner surface of the prosthesis. Brushes are tapered to provide effective cleaning of mucus and secretion build-up in the inner lumen of the voice prosthesis.



PB/50

## Provox

The Provox brush cleans the lumen and valve of the Provox voice prosthesis. Safety wings on the brush shaft help to prevent accidental aspiration of the Provox brush.



LAPCB0003

## Voiceline

A simple yet essential tool for cleaning leading brands of voice prosthesis. Specially constructed to help remove mucus and secretion build up on the inner surface of the prosthesis. Soft nylon bristles to help prevent damaging the inner silicone surface of the prosthesis.

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> IN4064	Blom Singer 16/20 Fr	Small 4mm - 8mm	3
<input type="checkbox"/> IN4065	Blom Singer 16/20 Fr	Medium 10mm - 16mm	3
<input type="checkbox"/> IN4066	Blom Singer 16/20 Fr	Large 18mm & above	3
<input type="checkbox"/> PB/50	Insight	-	10
<input type="checkbox"/> 7204	Provox	-	6
<input type="checkbox"/> 7225	Provox	X Large	6
<input type="checkbox"/> LAPCB0003	Voiceline	-	10

### Insight Tracheostomy Cannula Cleaning Swabs

For cleaning all types of PVC and silicone tracheostomy or laryngectomy tubes. The soft foam heads avoid damaging or scuffing the inside of the tube. Use with warm water to remove mucus residue in the tube.



**Check with your Head and Neck Nurse, ENT Nurse or Speech Therapist for suitability**



### Kapitex Trachi-Swab

Ideal for cleaning all PVC/PUC and silicone tracheostomy and laryngectomy tubes. The soft smooth heads of the swabs will not scratch or damage the internal tube surface.



Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/>	TCS/950 Insight Tracheostomy Cannula Cleaning Swabs	Small	10
<input type="checkbox"/>	TCS/900 Insight Tracheostomy Cannula Cleaning Swabs	Large	10
<input type="checkbox"/>	TR SWB 0150 Kapitex Trachi-Swab	Small (up to 5mm ID)	50
<input type="checkbox"/>	TR SWB 0250 Kapitex Trachi-Swab	Large (6mm to 8mm ID)	50
<input type="checkbox"/>	TR SWB 3050 Kapitex Trachi-Swab	X Large (over 8mm ID)	50
<input type="checkbox"/>	8250 Provox Swabs	Small	50
<input type="checkbox"/>	8251 Provox Swabs	Medium	50
<input type="checkbox"/>	8252 Provox Swabs	Large	50
<input type="checkbox"/>	8258 Provox Swabs	XL	50

**Optifahl Stoma Cleaning Wipes**



Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> 33200	Optifahl Stoma Cleaning Wipes	-	30
<input type="checkbox"/> 33260	Optifahl Stoma Cleaning Wipes	Tub	60

**Opticlear Adhesive Remover Wipes**




33500

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> PTA600	SpiroTect Baseplate Adhesive Remover Wipe	-	30
<input type="checkbox"/> 8012	Provox Adhesive Remover Wipes	-	50
<input type="checkbox"/> 33500	Opticlear Adhesive Remover Wipes	-	30

**Provox Silicone Glue**

(For use with the Provox Flexiderm Adhesive Baseplate)

Provox Silicone Glue is a liquid glue that can be used to improve the adhesion between the skin and the adhesive.



7720



	Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/>	7720	Provox Silicone Glue (for use with Flexiderm Adhesive Baseplate)	-	1

### Advadraw (T) Dressing



CR4416

**Check with your Head and Neck Nurse, ENT Nurse or Speech Therapist for suitability**

### Insight Advazorb Plus (T)

A hydrophilic foam tracheostomy dressing with incredible absorbency ability. Newly improved with a film backing added to ensure improved strength and absorbance for high/medium exuding wounds.



CR3900

### Permafoam Dressing 8x8cm Fenestrated

The hydrophilic foam dressing PermoFoam tracheostomy Y-cut for the supply of entry points such as tracheostomy tubes and special applications.

	Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/>	CR4416	Advadraw (T) Dressing	8cm x 8cm	20
<input type="checkbox"/>	CR4212	Insight Advazorb (T) Lite	8cm x 8cm	20
<input type="checkbox"/>	CR3900	Insight Advazorb Plus (T)	8cm x 8cm	20
<input type="checkbox"/>	23094	Metalline	8cm x 9cm	50
<input type="checkbox"/>	ACTIVHEALTRAC10X10CM	Activheal Tracheostomy 10x10cm	10cmx10cm	10
<input type="checkbox"/>	PERMAFOAM 8X8CM	Permafoam Dressing 8x8cm Fenestrated	8cmx8cm	10



### Sensotrach Alu Slit

30030



### Sensotrach Duo Slit

30609



	Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/>	30030	Sensotrach Alu Slit	8cm x 10cm	10
<input type="checkbox"/>	30609	Sensotrach Duo Slit	9cm x 9.8cm	10
<input type="checkbox"/>	30521	Sensotract Uno Slit	9cm x 9.8cm	10

### Insight Two-piece Latex Free

Two-piece design allows secure fitting for most sizes. Cough reflex section that allows patient movement and comfort while still holding the tube securely in position. Suitable for all types of tracheostomy tubes.



### Marpac Tube Holder with twill ties

One-piece tube holder made with comfortable foam laminate and extra long cotton twill.

### Marpac Collar with twill ties

Popular two-piece adjustable tube holder made with comfortable foam laminate and cotton twill/velcro fastener. Offers convenience of a velcro tie, but with soft, strong cotton twill sitting in the eyelets of the tube.

### Nex Tracheostomy Tube Holder

Laminated foam fabric for comfort and elasticity. Soft backing to prevent chafing and discomfort. One size fits all adults, suitable for all brands of tube.


### Optiflausch K 24-45cm

A wide, soft tube holder for improved comfort. The velcro fastener enables it to be easily and quickly adjusted to the required neck size.

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> TH100	Insight Two-piece Latex Free	Adult	10
<input type="checkbox"/> TH150	Insight Two-piece Latex Free	Child	10
<input type="checkbox"/> 106	Marpac Tube Holder with Twill Ties	Neo-Natal	25
<input type="checkbox"/> 107D	Marpac Tube Holder with Twill Ties	Paediatric	25
<input type="checkbox"/> 201	Marpac Collar with Twill Ties	Adult	10
<input type="checkbox"/> 203	Marpac Collar with Twill Ties	Neo-Natal/Paediatric	10
<input type="checkbox"/> HNTTHA	Nex Tracheostomy Tube Holder	-	10
<input type="checkbox"/> 32550	Optiflausch K	24cm - 45 cm	1

### SpiroTect Tube holder One-piece

Made from soft, skin-friendly material to avoid any discomfort or irritation to the neck. Easily adjustable to suit each patient with simple to use velcro closures.



PTA200

### SpiroTect Tube holder Two-piece

Two-piece for additional secure positioning.



PTA 201

### SpiroTect Tube holder Two-piece (Child)


Child version to comfortably fit children's necks.

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> PTA200	SpiroTect Tube Holder One-piece	-	10
<input type="checkbox"/> PTA201	SpiroTect Tube Holder Two-piece	Adult	10
<input type="checkbox"/> PTA301	SpiroTect Tube Holder Two-piece	Child	10
<input type="checkbox"/> 1652	Freeevent Neckband 2 piece	Small	10
<input type="checkbox"/> 1662	Freeevent Neckband 2 piece	Large	10
<input type="checkbox"/> 1651	Freeevent Neckband 1 piece	Small	10
<input type="checkbox"/> 1661	Freeevent Neckband 1 piece	Large	10



### Laryvox Simply Fix

Canula tube holders which assure the secure fixation of tracheostomy tubes attached variably to the neck area. Easily and quickly adjusted to the desired position



### Provox TubeHolder

Provox TubeHolder has been developed for use with Provox Life LaryTube and Provox Life LaryButton. It is suitable for users with sensitive skin around the stoma. The integrated clip connectors allow for optimal fit to the wings of the Provox Life LaryTube and LaryButton.



Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> 32680	Laryvox Simply Fix	-	40
<input type="checkbox"/> 32630	Neckfix Duo K	-	1
<input type="checkbox"/> 507800-000000	Rusch Neckband Paediatric	-	1
<input type="checkbox"/> 507900-000000	Rusch Neckband Adult	-	1
<input type="checkbox"/> 7668	Provox TubeHolder	-	1

### Provox LaryClip

This is a two-piece system that helps to optimise the airtight attachment of Provox LaryButton and Provox LaryTube to the stoma.

**Check with your Head and Neck Nurse, ENT Nurse or Speech Therapist for suitability**

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> 7669	Provox LaryClip	-	8



### Sofshield Laryngectomy Protector

Sofshield uses a unique toggle system that eliminates the need to tie cords. Pull cords for ease of use and comfortable tension for patient.



Check with your Head and Neck Nurse, ENT Nurse or Speech Therapist for suitability

AS3698



AS4051



AS/4056



### SpiroTect Protector Antibacterial

High porous foam increases airflow and reduces breathing resistance, helping to restore normal airway resistance. Excellent filtration, moistening and warming of air to protect the trachea.



PTA700



PTA702

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> AS3698	Sofshield Laryngectomy Protector White	Small	10
<input type="checkbox"/> AS3669	Sofshield Laryngectomy Protector White	Large	10
<input type="checkbox"/> AS3778	Sofshield Laryngectomy Protector Beige	Small	10
<input type="checkbox"/> AS3765	Sofshield Laryngectomy Protector Beige	Large	10
<input type="checkbox"/> AS3777	Sofshield Laryngectomy Protector Blue	Small	10
<input type="checkbox"/> AS3764	Sofshield Laryngectomy Protector Blue	Large	10
<input type="checkbox"/> AS4056	Softshield Laryngectomy Protector Light Blue	Medium	10
<input type="checkbox"/> AS4051	Sofshield Laryngectomy Protector Black	Large	10
<input type="checkbox"/> PTA700	SpiroTect Protector Antibacterial	Small	10
<input type="checkbox"/> PTA702	SpiroTect Protector Antibacterial	Large	10



## SpiroTect Protector


SpiroTect™ Protector comes in a range of colours and a choice of two sizes, large and small. They are specifically designed to protect Laryngectomy and Tracheostomy patients airway to give them the confidence to do everyday tasks.

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> PTA704	SpiroTect Protector Aqua	Small	10
<input type="checkbox"/> PTA708	SpiroTect Protector Aqua	Large	10
<input type="checkbox"/> PTA705	SpiroTect Protector Grey	Small	10
<input type="checkbox"/> PTA709	SpiroTect Protector Grey	Large	10
<input type="checkbox"/> PTA706	SpiroTect Protector Purple	Small	10
<input type="checkbox"/> PTA710	SpiroTect Protector Purple	Large	10
<input type="checkbox"/> PTA703	SpiroTect Protector White	Small	10
<input type="checkbox"/> PTA707	SpiroTect Protector White	Large	10

## Provox Protector

A thick, absorbing stoma protector for users with a large amount of secretions.

- Thick 3 layered protector
- Absorbing, soft and comfortable material
- Easy to close soft locking mechanism



Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> 7385	Provox Protector White	Small	10
<input type="checkbox"/> 7386	Provox Protector White	Large	10
<input type="checkbox"/> 7387	Provox Protector Slim White	Small	10
<input type="checkbox"/> 7389	Provox Protector Slim White	Large	10
<input type="checkbox"/> 7388	Provox Protector Slim Blue	Small	10
<input type="checkbox"/> 7390	Provox Protector Slim Blue	Large	10
<input type="checkbox"/> 7392	Provox Protector Air Blue	Small	10
<input type="checkbox"/> 7394	Provox Protector Air Blue	Large	10
<input type="checkbox"/> 7391	Provox Protector Air White	Small	10
<input type="checkbox"/> 7393	Provox Protector Air White	Large	10

**Kapitex Buchanan Lite**

Delivering proven filtration and heat moisture exchange benefits. Buchanan® Lite – combining essential airway protection and comfort!



LA BLT 2002



LA BLT 2003

**Kapitex Buchanan Protector**

Hydrolox 1 Foam is effective in removing dust and dirt particles in the air, as well as providing moisture and warmth to inhaled air. Pure cotton, sandwiching Hydrolox 1 Foam, allows a soft fabric feel against the neck.



LA BUP 0002



LA BUP 0001

**Kapitex Buchanan DeltaNex Protector**

The DeltaNex System offers all the benefits of the Buchanan protectors combined with touch to close fasteners for easier application.

**Hirst Protector**

A bib comprising of two layers of knitted net material which fastens at the rear of the neck with Velcro type fastener.

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> LA BLT 2001	Kapitex Buchanan Lite White	Small	10
<input type="checkbox"/> LA BLT 2002	Kapitex Buchanan Lite Blue	Small	10
<input type="checkbox"/> LA BLT 2003	Kapitex Buchanan Lite Beige	Small	10
<input type="checkbox"/> LA BLT 1001	Kapitex Buchanan Lite White	Large	10
<input type="checkbox"/> LA BLT 1002	Kapitex Buchanan Lite Blue	Large	10
<input type="checkbox"/> LA BLT 1003	Kapitex Buchanan Lite Beige	Large	10
<input type="checkbox"/> LA BUP 0002	Kapitex Buchanan Protector	Small	10
<input type="checkbox"/> LA BUP 0001	Kapitex Buchanan Protector	Large	10
<input type="checkbox"/> LA DNX 0001	Kapitex Buchanan DeltaNex Protector	Small	10
<input type="checkbox"/> HIRST BIB	Hirst Protector	-	10



## Laryngofoam

Laryngofoam filters offer the laryngectomee lightweight, comfortable and unobtrusive protection. The filters warm, moisten and filter the incoming air and can be applied in seconds.



**Check with your Head and Neck Nurse, ENT Nurse or Speech Therapist for suitability**

LALFF0023

## Purifoam Laryngectomy Filter

Purifoam uses a high quality medical grade foam that protects the stoma from dirt and dust entering the airway, while warming and humidifying the incoming air.

## Provox Tracheofix

Provox Tracheofix is a single-use 55 x 60 mm (small) foam protector that provides aesthetic coverage of tracheostoma. Helps to prevent dust and dirt from entering the stoma.



1433H

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> LALFF0003	Laryngofoam Flesh	5.1cm x 6.2cm	30
<input type="checkbox"/> LALFF0023	Laryngofoam Flesh	6.3cm x 6.5cm	30
<input type="checkbox"/> CR3853	Purifoam Laryngectomy Filter	6.5cm x 6.5cm	30
<input type="checkbox"/> 1429H	Provox Tracheofix (1 strip) Beige	5.5cm x 6cm	10
<input type="checkbox"/> 1432H	Provox Tracheofix (2 strip) Beige	5.5cm x 6cm	10
<input type="checkbox"/> 1430H	Provox Tracheofix (1 strip) Beige	7cm x 7cm	10
<input type="checkbox"/> 1433H	Provox Tracheofix (2 strip) Beige	7cm x 7cm	10
<input type="checkbox"/> 1427	Provox Tracheofix (1 strip) Ivory	5.5cm x 6cm	10
<input type="checkbox"/> 1428	Provox Tracheofix (1 strip) Ivory	7cm x 7cm	10

## Cascade Shower Protector

Cascade is a soft, waterproof PVC protector that provides a barrier to prevent shower water from entering a stoma thereby improving the patient's confidence to shower routinely.



AS/3835



Check with your Head and Neck Nurse, ENT Nurse or Speech Therapist for suitability

## Spiroprotect Shower Protector

Enables a laryngectomee to shower without water entering the trachea. It is small and simple to use, taking the hassle out of showering with a Laryngectomy stoma.



PTA101

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> AS/3835	Cascade Shower Protector	-	2
<input type="checkbox"/> PTA101	Spiroprotect Shower Protector	-	2
<input type="checkbox"/> LAPPB0001	Phoniq DigiClose Shower Shield	-	1
<input type="checkbox"/> LAPPB0002	Phoniq Low Profile Shower Shield	-	2
<input type="checkbox"/> LAPPB0003	Phoniq Universal Shower Shield	-	2
<input type="checkbox"/> 8308	Provox Life Shower	-	1
<input type="checkbox"/> 8016	Provox Luna Shower Aid	-	1
<input type="checkbox"/> 7260	Provox Shower Aid	-	1



## Barrier Cream

Skin barriers protect the skin from moisture associated skin damage and reduce the risk of maceration and skin stripping from adhesive devices.



Check with your Head and Neck Nurse, ENT Nurse or Speech Therapist for suitability

### Cavilon Barrier Cream and Film



3392GS



3392G

### Sorbaderm Barrier Cream



3028

### SpiroTect Baseplate Barrier Wipe

Leaves a protective and quick drying layer on the skin. Helps to reduce skin irritation from baseplates. The wipe contains a special non-sting formula and does not contain alcohol.




PTA500

	Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/>	3392G	Cavilon Barrier Cream	92g	1
<input type="checkbox"/>	3392GS	Cavilon Barrier Cream	2g	20
<input type="checkbox"/>	3343P	Cavilon Barrier Film	1ml	5
<input type="checkbox"/>	3028	Sorbaderm Barrier Cream	92g	1
<input type="checkbox"/>	33600	Optiguard Protective Dressing Wipes	-	30
<input type="checkbox"/>	PTA500	SpiroTect Baseplate Barrier Wipe	-	30
<input type="checkbox"/>	7244	Provox Cleaning Towel	-	10
<input type="checkbox"/>	8011	Provox Skin Barrier Wipes	-	50

### Alvita Saline Pods



### Normasol Irrigation Solution



Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> ALVITA SALINE PODS	Alvita Sterile Saline Pods	20ml	25
<input type="checkbox"/> IRR24N	Irripod	20ml	25
<input type="checkbox"/> NOR206B	Normasol	25ml	25

### Alvita Gauze Non Sterile 10cm x 10cm 8 ply Gauze Swabs



GAUZE810

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> GAUZE810	Alvita Gauze Non Sterile 10cm x 10cm 8 ply swabs	-	100

### Alvita Sterile Dressing Pack Spec 10

Spec 10 pack contains a gauze and cotton tissue pad, gauze swabs, cotton balls, paper towel and an inner wrapping to be used as a sterile working field.



STEV10 ALVITA

### Alvita Sterile Dressing Pack Spec 35

Spec 35 pack contains a non woven fabric covered dressing pad, nonwoven fabric swabs, cotton wool balls, paper towel and an inner wrapping to be used as a sterile working field.



STEV35 ALVITA

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> STEV10 ALVITA	Sterile Dressing Pack Spec 10	-	12
<input type="checkbox"/> STEV35 ALVITA	Sterile Dressing Pack Spec 35	-	12

### Optilube Lubricating Jelly

Optilube lubricating jelly is a water based, water soluble gel that is safe to use and aids insertion of tubes and studs.



1120



1122

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> 1120	Optilube (sachet)	5g	150
<input type="checkbox"/> 1121	Optilube (tube)	42g	1
<input type="checkbox"/> 1122	Optilube (tube)	82g	1

### Provox Adhesive Strips

Provox Adhesive Strips can be used to cover the edges of a baseplate to help prevent it from loosening from the skin.

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> 8015	Provox Adhesive Strips	-	10

### Insil Soft Silicone Tape

Insil is a versatile tape backed with a silicone adhesive. When the tape is removed it is gentle on the skin, however it can be re-applied many times without losing any adhesive qualities.



INSIL400

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> INSIL200	Insil	2cm x 3m	1
<input type="checkbox"/> INSIL400	Insil	4cm x 1.5m	1

LATNG2008



LATNG4008



### Kapi-Gel Spacers

Kapi-Gel Spacers are ideal for difficult stoma anatomy and help to fill the stoma area, allowing the comfortable use of Laryngectomy tubes, stoma studs and HME/filter holders.

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> LATNG2008	Kapi-Gel Spacers Thin	3mm	10
<input type="checkbox"/> LATNG2012	Kapi-Gel Spacers Thin	3mm	10
<input type="checkbox"/> LATNG4008	Kapi-Gel Spacers Thick	5mm	10
<input type="checkbox"/> LATNG4012	Kapi-Gel Spacers Thick	5mm	10



### TheraBite System Jaw Rehabilitation

TheraBite Jaw Motion Rehabilitation System can help to reduce muscle pain by reducing joint inflammation and preventing reduced muscle strength during (chemo) radiation in the head and neck area.

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> TH001	TheraBite Jaw Motion Rehabilitation System	Adult	1
<input type="checkbox"/> TH002	TheraBite Jaw Motion Rehabilitation System	Paediatric	1
<input type="checkbox"/> PA001	TheraBite Bite Pad	Regular	4
<input type="checkbox"/> PA002	TheraBite Bite Pad	Endentulous	4

### Gelclair Oral Treatment

Gelclair oral spray and sachets.

Product Code	Item Description	Size/Shape	Pack Size
<input type="checkbox"/> GelX001	Oral Spray for prevention and treatment of oral Mucositis		1
<input type="checkbox"/> Gelclair	Gelclair Sachets	21 x 15ml	21







## How To Use The NWOS Homecare Delivery Service



Request and collect your prescription from your Doctor/specialist. Check your details and requirements are correct. Sign the back of the prescription. We can do this for you with your permission.



Send us your order by: email [homecareorders@nwossurgical.co.uk](mailto:homecareorders@nwossurgical.co.uk), SMS **07464 986 747** or Freephone **0800 243 103** (between **8.30am** and **5.00pm**)



Post your prescription to us in the pre-paid envelope provided. Please send a copy of your exemption certificate if you are using us for the first time.



We will deliver your order direct to your door.



Remember:

- Order your prescription in good time if you are on holiday or going away
- If you have any specific delivery instructions e.g. if you would like us to leave your order in a garage, shed or an alternative address, then please let us know
- If you telephone us in advance, we can deliver your order to an alternative address in the UK

### Contact us:



Email [homecareorders@nwossurgical.co.uk](mailto:homecareorders@nwossurgical.co.uk)



Freephone **0800 243 103** 8.30am to 5.00pm



SMS text number **07464 986 747**



[www.homecarenwos.co.uk](http://www.homecarenwos.co.uk)



**Homecare**<sup>®</sup>  
from NWOS<sup>®</sup>



Free delivery



See above for full details



**Leigh Service Centre**  
Green Fold Way,  
Leigh Commerce Park,  
Leigh, WN7 3XJ

